

PORTABLE OXYGEN CONCENTRATORS

FAA APPROVED OCTOBER 31, 2012



RENTALS, SALES, SUPPORT
AEROMEDIC.COM
1-800-346-3556

PICTURES ARE NOT TO SCALE. CHECK YOUR AIRLINES WEBSITE FOR RULES ABOUT FLYING WITH PORTABLE OXYGEN CONCENTRATORS.



AirSep Focus™
 Weight: 1.75 lbs
 Dims: 6.2" x 4.6" x 2.4"



AirSep FreeStyle™
 Weight: 4.4 lbs
 Dims: 8.6" x 6.1" x 3.6"



AirSep FreeStyle 5™
 Weight: 6.2 lbs
 Dims: 10.7" x 6.6" x 4.4"



Invacare XPO2™
 Weight: 6 lbs
 Dims: 10" x 7" x 4"



Inogen One™
 Weight: 9.7 lbs
 Dims: 11.62" x 6" x 12.39"



Inogen G2™
 Weight: 7 lbs
 Dims: 11.61" x 5.98" x 10.75"



AirSep LifeStyle™
 Weight: 9.75 lbs
 Dims: 5.5" x 7.25" x 16.31"



Precision Medical EasyPulse™
 Weight: 6.8 lbs
 Dims: 10.1" x 6.5" x 4.5"



Inova Labs LifeChoice™
 Weight: 4.9 lbs
 Dims: 9.5" x 7.5" x 3.125"



Inova Labs ActiveOx™
 Weight: 4.83 lbs
 Dims: 9.05" x 7.875" x 4.38"



Inogen G3™
 Weight: 4.8 lbs
 Dims: 8.75" x 3.0" x 7.25"



SeQual Eclipse®
 Weight: 17.9 lbs
 Dims: 19.3" x 12.3" x 7.1"



Philips Respironics EverGo™
 Weight: 8.5 lbs
 Dims: 12" x 6" x 8.5"



Oxus™
 Weight: 9.8 lbs
 Dims: 11.6" x 7.4" x 4.6"



DeVilbiss Healthcare iGo™
 Weight: 19 lbs
 Dims: 15" x 11" x 8"



Invacare Solo2™
 Weight: 19.9 lbs
 Dims: 16.5" x 11" x 8"



SeQual SAROS®
 Weight: 12.25 lbs
 Dims: 26.8" x 4.375"



Philips Respironics SimplyGo™
 Weight: 10 lbs
 Dims: 11.5" x 10" x 6"



OxLife Independence™
 Weight: 16.7 lbs
 Dims: 20.29" x 10.85" x 9.45"

CONSULT YOUR DOCTOR ABOUT FLYING WITH OXYGEN. NEED MORE CARDS LIKE THESE? EMAIL INFO@AEROMEDIC.COM

PORTABLE OXYGEN CONCENTRATORS

BATTERY DURATION

FAA APPROVED

OCTOBER 31, 2012

AEROMEDIC.COM

1

2

3

Setting

4

5

6

AirSep Focus	--	1.5 hrs	--	--	--	--
AirSep Battery Belt +add	--	4.0 Hrs	--	--	--	--
AirSep FreeStyle	3.5 hrs	2.5 hrs	2.0 hrs	--	--	--
AirSep Battery Belt +add	6.5 hrs	3.5 hrs	3.0 hrs	--	--	--
AirSep FreeStyle 5	3.0 hrs	2.0 hrs	1.5 hrs	1.25 hrs	1.0 hrs	--
AirSep Battery Belt +add	4.0 hrs	2.5 hrs	2.0 hrs	1.75 hrs	1.5 hrs	--
AirSep LifeStyle	50 min	50 min	50 min	50 min	50 min	--
Inogen One	3.0 hrs	3.0 hrs	2.5 hrs	2.5 hrs	2.0 hrs	--
Inogen One G2 12 cell	4.0 hrs	3.5 hrs	2.7 hrs	2.2 hrs	2.0 hrs	--
G2 24 cell Extended Battery	8.0 hrs	7.0 hrs	5.4 hrs	4.4 hrs	4.0 hrs	--
Inogen One G3 8 cell	4.5 hrs	4.0 hrs	3.0 hrs	2.0 hrs	--	--
G3 16 cell Extended Battery	9.0 hrs	8.0 hrs	6.0 hrs	4.0 hrs	--	--
Inova Labs LifeChoice	2.0 hrs	2.0 hrs	2.0 hrs	--	--	--
Lifechoice External Battery +add	3.0 hrs	3.0 hrs	3.0 hrs	--	--	--
Inova Labs LifeChoice Activox	12.0 hrs	6.0 hrs	5.0 hrs	--	--	--
Activeox External Battery +add	3.0 hrs	3.0 hrs	2.0 hrs	--	--	--
Invacare XPO2	3.5 hrs	2.5 hrs	2.0 hrs	1.5 hrs	1.0 hrs	--
XPO2 External Battery +add	3.5 hrs	2.5 hrs	2.0 hrs	1.5 hrs	1.0 hrs	--
Oxus RS-00400, Delphi RS-00400	3.45 hrs	3.0 hrs	2.5 hrs	2.0 hrs	1.75 hrs	--
Precision Medical EasyPulse	4.5 hrs	3.2 hrs	2.5 hrs	1.8 hrs	1.5 hrs	--
EasyPulse External Battery +add	4.5 hrs	3.2 hrs	2.5 hrs	1.8 hrs	1.5 hrs	--
Philips Respironics EverGo	6.0 hrs	4.0 hrs	3.0 hrs	2.0 hrs	2.0 hrs	2.0 hrs
EverGo with 2nd Battery	12.0 hrs	8.0 hrs	6.0 hrs	4.0 hrs	4.0 hrs	4.0 hrs
DeVilbiss Healthcare iGo	5.4 hrs	4.7 hrs	4.0 hrs	3.5 hrs	3.2 hrs	3.0 hrs
Continuous Flow	4.4 hrs	2.4 hrs	1.6 hrs	--	--	--
Invacare Solo2*	4.5 hrs	3.5 hrs	3.0 hrs	2.5 hrs	2.5 hrs	--
Continuous Flow	3.5 hrs	2.5 hrs	1.5 hrs	--	--	--
Oxlife Independence*	2.8 hrs	2.8 hrs	1.4hrs	1.4 hrs	53 min	--
With 2nd Battery	5.6 hrs	5.6 hrs	2.8 hrs	2.8 hrs	1.8 hrs	--
Continuous Flow Battery / w/ 2nd Battery	2.4/4.8 hrs	1.4/2.8 hrs	53min/1.8hr	--	--	--
Philips Respironics SimplyGo*	3.7 hrs	3.5 hrs	3.1 hrs	2.7 hrs	2.3 hrs	1.9 hrs
Continuous Flow	2.3 hrs	0.7 hrs	--	--	--	--
SeQual Eclipse*	4.4 hrs	3.6 hrs	3.0 hrs	2.6 hrs	2.3 hrs	2.1 hrs
Continuous Flow	3.7 hrs	2.0 hrs	1.3 hrs	--	--	--
SeQual SAROS	1.2 hrs	1.1 hrs	53 min	47 min	46 min	37 min
Continuous Flow	1.1 hrs	45 min	30 min	--	--	--

All times based on NEW batteries. Each airline has there own rules about flying with portable oxygen concentrators *Solo2 (.5 lpm = 4.5 hrs), *Independence (.5 lpm = 2.4 hrs) *SimplyGo (.5 lpm = 3.1 hrs), *Eclipse (.5 lpm = 4.4 hrs) A doctors prescription is required for all portable oxygen concentrators. Each airline has there own specific Physician Statement to be completed, doctor signed, and carried with you on the plane. Talk to your doctor about Pulse or On-Demand Delivery of Oxygen. Advanced Aeromedical is On-Call 24 hrs for support. Rentals / Sales / Support 1-800-346-3556 info@aeromedical.com



O₂